


# Eliminates the Dead Spot.



## ABILITY DYNAMICS™

### RUSH™ Foot Single-Step Analysis



### CATEGORY CHART

Weight (lbs)	90 – 116	117 – 146	147 – 180	181 – 218	219 – 261	262 – 310	311 – 360
Low Activity	1	2	3	4	5	6	7
Moderate Activity	1	2	3	4	5	6	7
High Activity	2	3	4	5	6	7	8
High Impact	2-3	3-4	4-5	5-6	6-7	7-8	8-9

NOTE: Increasing category numbers indicate increased keel response.

Call toll-free at 855-450-7300  
or visit us at [abilitydynamics.com](http://abilitydynamics.com).

Proudly handcrafted in the U.S.A.

