

MEET THE NEW **KID** ON THE BLOCK

ALL THE ABILITIES OF THE RUSH87™ IN JUST UNDER 2"!



- SMOOTH ROLL-THROUGH
- DYNAMIC TOE-OFF
- MOST DYNAMIC LOW PROFILE FOOT AVAILABLE
- THE ONLY LOW PROFILE 81/86 FOOT ON MARKET

PROFILE INFORMATION

Weight Limit	125lbs/57kg	Heel Rise	3/8"/10mm
Size	19cm – 21cm	Suggested Codes	L-5981/5986
Foot Shell	Light or Dark	Warranty	18 months
Product Weight	276g	Build Height	1.75"*

* Based on 21cm, Category 4P.

CATEGORY CHART

Weight (in lbs)	35-56	57-78	79-101	102-125
Moderate Activity	1P	2P	3P	4P
High Activity	2P	3P	4P	5P
High Impact	2P	3P	4P	5P

Call us toll-free at 855.450.7300 or visit us at rushfoot.com

f t YouTube in p #notjustanothercarboncopykid



Proudly handcrafted in the U.S.A.

RUSHTM
Kid

