

RUSH™ ROGUE



INCREASED PERFORMANCE. INCREASED COMFORT. ADDED TORSION.

RUSH™ ROGUE provides the most realistic and dynamic foot and ankle motion available. The Vertical Loading Pylon provides comfortable relief from high impact loading allowing the user to push the foot to the limits. The components are made of advanced fiberglass composite which is three times more flexible than most conventional prosthetic feet. The roll through characteristics of the foot provide exceptional energy return with no “dead spot” allowing users to push themselves to the limits on the court or in aggressive terrain.

- **Vertical Loading Pylon** for vertical shock relief and impact absorption
- **Smooth roll through** for high energy return and comfortable in any situation
- **Biomimetic ankle action** to simulate the human foot – no one will know you wear a prosthetic
- **Extreme reliability** for extended use in any high impact application a human can survive
- **Rotational torsion relief** to ease the impact to knee, hip and back



PROFILE INFORMATION

PROFILE INFORMATION			BUILD HEIGHT		
Weight Limit	360lbs./150kg	Heel Rise	3/8”/10mm	Size	Build Height*
Size	22cm - 29cm	†Suggested Codes	L5987, L5984	22 - 24cm	6.5”
Foot Shell	Light or Dark	Warranty	36 Months	25 - 27cm	6.8”
Product Weight	865g/2.18lbs. (26cm)			28 - 29cm	7.3”

† The listing of this code should not be considered a guarantee for coverage or a payment. * Build height includes footshell.

CATEGORY CHART

Weight (lbs)	90 - 116	117 - 146	147 - 180	181 - 218	219 - 261	262 - 310	311 - 360
Low Activity	1	2	3	4	5	6	7
Moderate Activity	1	2	3	4	5	6	7
High Activity	2	3	4	5	6	7	8
High Impact	2-3	3-4	4-5	5-6	6-7	7-8	8-9

NOTE: Increasing category numbers indicate increased keel response.

SMOOTH · DYNAMIC · RELIABLE · VIRTUALLY INDESTRUCTIBLE