



# RUSH Chopart®

## VIRTUALLY INDESTRUCTIBLE

The **RUSH Chopart®** is specifically designed for people with Chopart amputations as well as lisfranc, pirogoff, boyd and partial foot prostheses. A dynamic design, providing increased reliability and eliminating the breakage commonly seen with carbon Chopart plates. The unique glass composite is engineered to be more flexible and durable than any other Chopart available today.



### PROFILE INFORMATION

<b>Weight Limit</b>	360lbs./163kg	<b>Heel Rise</b>	3/8"/10mm	<b>Build Height*</b>
<b>Size</b>	22cm - 30cm	<b>Warranty</b>	24 Months	.285" + 3/8"* (* based on 27cm Cat 4)
<b>Foot Shell</b>	Light or Dark			
<b>Product Weight</b>	137g* (* based on 27cm Cat 4)			

† The listing of this code should not be considered a guarantee for coverage or a payment. \* Build height includes footshell.

### CATEGORY CHART

Weight	lbs	90 - 116	117 - 146	147 - 180	181 - 218	219 - 261	262 - 310	311 - 360
	kg	40-52	53-66	67-81	82-98	99-118	119-140	141-163
<b>Low Activity</b>		1	2	3	4	5	6	7
<b>Moderate Activity</b>		1	2	3	4	5	6	7
<b>High Activity</b>		2	3	4	5	6	7	8
<b>High Impact</b>		2-3	3-4	4-5	5-6	6-7	7-8	8-9

NOTE: Increasing category numbers indicate increased keel response.

Call us toll-free at 855.450.7300 or visit us at [proteorusa.com](http://proteorusa.com)

f t You Tube in p #HumanFirst



HumanFirst