

SIMPLISTIC. EFFECTIVE. EFFICIENT.

Finally, an unprecedented elevated vacuum system that seamlessly integrates effective vacuum with the same rugged simplicity of all the **RUSH Foot®** designs. The **RUSH Foot EVAQ8 Collection** is the easy solution to providing maximum vacuum efficiency in just a few quick steps. No batteries or complicated, mechanical pump required.

Unlike the competition, upwards of 22inHg can be achieved with this surprisingly simple, integrated design that stabilizes volume, increases linkage, and heightens proprioception through elevated vacuum suspension. The **RUSH Foot EVAQ8 foot Collection** is all-terrain, fully-submersible (fresh or salt water), free of bulky components, and requires minimal maintenance. The **RUSH Foot EVAQ8 Collection** is available in the **RUSH LoPro**, **HiPro** and **ROGUE** models providing the perfect option for the majority of individuals seeking elevated vacuum suspension.

Go anywhere. Do anything. The RUSH Foot EVAQ8 Collection elevates your life.

RUSH FOOT® EVAQ8 COLLECTION



RUSH ROGUE® EVAQ8 PROFILE INFORMATION

Weight Limit	360lbs./163kg
Size	22cm - 29cm
Foot Shell	Light or Dark
Product Weight	875g/1.93lbs. (26cm)

BUILD HEIGHT

Size	Build Height*
22 - 24cm	6.5"
25 - 27cm	6.8"
28 - 29cm	7.3"

RUSH HiPro® EVAQ8 PROFILE INFORMATION

Weight Limit	360lbs./163kg
Size	22cm - 29cm
Foot Shell	Light or Dark
Product Weight	770g/1.70lbs. (26cm)

BUILD HEIGHT

Size	Build Height*
22 - 24cm	6.5"
25 - 27cm	6.8"
28 - 29cm	7.3"

RUSH LoPro® EVAQ8 PROFILE INFORMATION

Weight Limit	300lbs./136kg
Size	22cm - 29cm
Foot Shell	Light or Dark
Product Weight	649g/1.43lbs. (26cm)

BUILD HEIGHT

Size	Build Height*
22 - 24cm	4.5"
25 - 27cm	4.5"
28 - 29cm	4.5"

† The listing of this code should not be considered a guarantee for coverage or a payment. * Build height includes footshell.

CATEGORY CHART

Weight (lbs)	90 - 116	117 - 146	147 - 180	181 - 218	219 - 261	262 - 310	262 - 300 LoPro Limit	311 - 360
Low Activity	1	2	3	4	5	6	7	7
Moderate Activity	1	2	3	4	5	6	7	7
High Activity	2	3	4	5	6	7	8	8
High Impact	2-3	3-4	4-5	5-6	6-7	7-8	8-9	8-9

NOTE: Increasing category numbers indicate increased keel response.

Call us toll-free at 855.450.7300 or visit us at proteorusa.com

f t You Tube in p #HumanFirst



HumanFirst