

SMOOTH · FLEXIBLE · RELIABLE · VIRTUALLY INDESTRUCTIBLE

Trenten Merrill, BK, RUSH HiPro®  
2016 Team USA Paralympic Athlete

# RUSH HiPro®

## VIRTUALLY INDESTRUCTIBLE

RUSH HiPro® provides the most realistic and responsive foot and ankle motion available. Made from an advanced engineered glass composite, which is three times more flexible than most conventional prosthetic feet. The energy return of this foot enables amputees to easily traverse the most aggressive terrain while experiencing zero dead spots.

**“THE IMPACT ABSORPTION IS UNREAL,  
IT WITHSTANDS EVERYTHING I CAN  
THROW AT IT.”**

– Jay Nero, BK, RUSH HiPro®



### PROFILE INFORMATION

<b>Weight Limit</b>	360lbs./163kg	<b>Heel Rise</b>	3/8"/10mm	<b>Size</b>	<b>Build Height*</b>
<b>Size</b>	22cm - 29cm	<b>Warranty</b>	36 Months	22 - 24cm	6.5"
<b>Foot Shell</b>	Light or Dark			25 - 27cm	6.8"
<b>Product Weight</b>	760g/1.67lbs. (26cm)			28 - 29cm	7.3"

† The listing of this code should not be considered a guarantee for coverage or a payment. \* Build height includes footshell.

### CATEGORY CHART

Weight	90 - 116	117 - 146	147 - 180	181 - 218	219 - 261	262 - 310	311 - 360
lbs	40-52	53-66	67-81	82-98	99-118	119-140	141-163
kg							
<b>Low Activity</b>	1	2	3	4	5	6	7
<b>Moderate Activity</b>	1	2	3	4	5	6	7
<b>High Activity</b>	2	3	4	5	6	7	8
<b>High Impact</b>	2-3	3-4	4-5	5-6	6-7	7-8	8-9

NOTE: Increasing category numbers indicate increased keel response.

Call us toll-free at 855.450.7300 or visit us at [proteorusa.com](http://proteorusa.com)

f t You Tube in p #HumanFirst



Human First