

Introducing the ultra low profile RUSH Foot® now in sizes 22-30cm!

All the abilities of the RUSH HiPro® in just over 2"!



- SMOOTH ROLL-THROUGH
- DYNAMIC TOE-OFF
- MOST DURABLE LOW-PROFILE FOOT AVAILABLE

PROFILE INFORMATION		BUILD HEIGHT			
Weight Limit	360lbs/163kg	Heel Rise	3/8"/10mm	Size	Build Height*
Size	22cm – 30cm	†Suggested Codes	L-5981/5986	22-24cm	2.1"
Foot Shell	Light or Dark	Warranty	36 months	25-27cm	2.2"
Product Weight	950g/2.1lbs			28-30cm	2.3"

* The listing of this code should not be considered a guarantee for coverage or a payment. * Build height includes footshell.

CATEGORY CHART							
Weight (in lbs)	90-116	117-146	147-180	181-218	219-261	262-310	311-360
Low Activity	1	2	3	4	5	6	7
Moderate Activity	1	2	3	4	5	6	7
High Activity	2	3	4	5	6	7	8
High Impact	2-3	3-4	4-5	5-6	6-7	7-8	8-9

*This information is preliminary and is subject to change at time of release.



Call us toll-free at 855.450.7300 or visit us at proteorusa.com

f t YouTube in p i #HumanFirst



Human First