



# RUSH™

## RAMPAGE

## Unbridled Capabilities.

The latest addition to the family, the **RUSH Foot RAMPAGE**, has arrived. We've taken the smooth roll-through, increased comfort, and all-terrain capability of its predecessors and created a new, more cosmesis-friendly option, with a lower overall product weight. The **RUSH Foot RAMPAGE** is the perfect solution for individuals seeking a blend of both cosmetic, and high-performance capabilities found only in the RUSH Foot Collection. Are you ready to **RAMPAGE?**

- Cosmesis-ready to meet your aesthetic preferences.
- Lower overall product weight.
- Smooth roll through with unparalleled energy return.
- Biometric ankle motion simulates natural gait, increasing patient comfort.
- Extreme reliability for any use in any environment.
- Manufacturer recommended coding based on AOPA testing guidelines L5987 and L5986

PROFILE INFORMATION		PRODUCT WEIGHT*		BUILD HEIGHT**	
Weight Limit	365lbs./165kg	Foot Only	1.32lbs./600g	Size	Build Height
Size	22cm - 29cm	With Foot Shell and Sock	1.97lbs./895g	22 - 24cm	6.5"
Foot Shell	Light or Dark			25 - 27cm	6.9"
Heel Rise	3/8"/10mm			28 - 29cm	7.3"
Warranty	36 Months	Lcodes L5987 and L5986			

\*based on size 27 Cat 4

\*\* Build height includes footshell.

### CATEGORY CHART

Weight	90 - 116	117 - 146	147 - 180	181 - 218	219 - 261	262 - 310	311 - 365
lbs	40-52	53-66	67-81	82-98	99-118	119-140	141-165
kg							
Low Activity	1	2	3	4	5	6	7
Moderate Activity	1	2	3	4	5	6	7
High Activity	2	3	4	5	6	7	8
High Impact	2-3	3-4	4-5	5-6	6-7	7-8	8-9

NOTE: Increasing category numbers indicate increased keel response.

Call us toll-free at 855.450.7300 or visit us at [proteorusa.com](http://proteorusa.com)

f t YouTube in p #HumanFirst



HumanFirst