

# RUSH™

## RAMPAGE



## Unbridled Capabilities.

The latest addition to the family, the **RUSH Foot RAMPAGE**, has arrived. We've taken the smooth roll-through, increased comfort, and all-terrain capability of its predecessors and created a new, more cosmesis-friendly option, with a lower overall product weight. The **RUSH Foot RAMPAGE** is the perfect solution for individuals seeking a blend of both cosmetic, and high-performance capabilities found only in the RUSH Foot Collection. Are you ready to **RAMPAGE?**

- Cosmesis-ready to meet your aesthetic preferences.
- Lower overall product weight.
- Smooth roll through with unparalleled energy return.
- Biometric ankle motion simulates natural gait, increasing patient comfort.
- Extreme reliability for any use in any environment.
- Manufacturer recommended coding based on AOPA testing guidelines L5987 and L5986

PROFILE INFORMATION		PRODUCT WEIGHT*		BUILD HEIGHT**		CATEGORY CHART								
Weight Limit	360lbs./168kg	Standard	1.29lbs./583g	Size	Build Height	Weight	lbs	90 - 116	117 - 146	147 - 180	181 - 218	219 - 261	262 - 310	311 - 365
Size	22cm - 29cm			22 - 24cm	6.5"	kg	40-52	53-66	67-81	82-98	99-118	119-140	141-165	
Foot Shell	Light or Dark	EVAQ8	1.31lbs./596g	25 - 27cm	6.9"	Low Activity		1	2	3	4	5	6	7
Heel Rise	3/8"/10mm	H2O	1.49lbs./674g	28 - 29cm	7.3"	Moderate Activity		1	2	3	4	5	6	7
Warranty	36 Months			Lcodes L5987 and L5986		High Activity		2	3	4	5	6	7	8
						High Impact		2-3	3-4	4-5	5-6	6-7	7-8	8-9

\*Based on size 26 Cat 4, foot only

\*\* Build height includes footshell.

NOTE: Increasing category numbers indicate increased keel response.

**NOW AVAILABLE IN EVAQ8 AND H2O MODELS!**

Call us toll-free at 855.450.7300 or visit us at [proteorusa.com](http://proteorusa.com)

f t YouTube in p i #HumanFirst



Human First