











# Which RUSH Foot® is Right for You?

Call 855.450.7300 or visit us at [rushfoot.com](http://rushfoot.com)

	DESCRIPTION	WEIGHT LIMIT	SIZES	IMPACT LEVELS	CLEARANCE	WARRANTY	FOOT SHELL	SUGGESTED CODING	PDAC APPROVED CODING
 <b>RUSH ROGUE®</b>	Vertical Loading Pylon Smooth roll through Biomimetic ankle action Extreme reliability Rotational torsion relief Waterproof—fresh or salt water	360 lbs	22cm-29cm	low, moderate, high, high impact	6.5"-7.3"	36 months	Light or Dark	L5987 + L5984 + L5986	L5987 + L5984
 <b>RUSH HiPro®</b>	Highly-dynamic Durable Smooth roll-over High-energy return All-terrain Waterproof—fresh or salt water Reduced socket issues	360 lbs	22cm-29cm	low, moderate, high, high impact	6.5"-7.3"	36 months	Light or Dark	L5987 + L5986	L5981
 <b>RUSH LoPro®</b>	Low-profile build Superior ground compliance Smooth roll-over High-energy return Cosmesis-friendly All-terrain Waterproof—fresh or salt water Reduced socket issues	360 lbs	22cm-29cm	low, moderate, high, high impact	4.5"	36 months	Light or Dark	L5981 + L5986	L5981
 <b>RUSH ROVER®</b>	Highly-dynamic Extremely flexible Dynamic toe-off Smooth roll-through Linear Heel Response Virtually indestructible Waterproof—fresh or salt water Reduced socket, knee and hip pain Cosmesis-friendly	360 lbs	22cm-30cm	low, moderate, high, high impact	2.1"-2.3"	36 months	Light or Dark	L5981 + L5986	L5981
 <b>RUSH Kid®</b>	Highly-dynamic Extremely flexible Dynamic toe-off Smooth roll-through Virtually indestructible Waterproof—fresh or salt water Reduced socket, knee and hip pain Cosmesis-friendly	125 lbs	19cm-22cm	moderate, high, high activity	1.75"	18 months	Light or Dark	L5981 + L5986	N/A
 <b>RUSH Chopart®</b>	Highly-dynamic Durable Smooth roll-over Waterproof—fresh or salt water	360 lbs	22cm-30cm	moderate, high, high activity	17mm (based on 27cm Cat 4)	24 months	Light or Dark	L5976	L5976
 <b>RUSH H2O</b>	All the same great features as our standard RUSH feet, but intended for water or barefoot activities only. No foot shell.	360 lbs	22cm-29cm	low, moderate, high, high impact	Rogue: 6.5"-7.3" HiPro: 6.5"-7.3" LoPro: 4.5"	36 months	No foot shell required	Same as the standard ROGUE, HiPro and LoPro	N/A
 <b>RUSH EVAQ8</b>	All the same great features as our standard RUSH feet plus integrated vacuum system.	360 lbs	22cm-29cm	low, moderate, high, high impact	Rogue: 6.5"-7.3" HiPro: 6.5"-7.3" LoPro: 4.5"	36 months	Light or Dark	Same as the standard ROGUE, HiPro and LoPro + L5781	N/A