

# RUSH<sup>FOOT</sup><sup>®</sup>

## EVAQ8

### COLLECTION

## AT-A-GLANCE VACUUM FOOT COMPARISON

	Requires Charging	Maximum Vacuum	Sleeve Requirement	Waterproof	Mounting Installation	Field Serviceable	Filter Serviceable
 <p><b>ROGUE EVAQ8</b> Manufacturer: Ability Dynamics</p>	NO	22 inHg	No sleeve required if used with seal-in liner	YES	INTEGRATED	YES	YES
<a href="http://rushfoot.com">rushfoot.com</a>							
 <p><b>HiPro EVAQ8</b> Manufacturer: Ability Dynamics</p>	NO	22 inHg	No sleeve required if used with seal-in liner	YES	INTEGRATED	YES	YES
<a href="http://rushfoot.com">rushfoot.com</a>							
 <p><b>LoPro EVAQ8</b> Manufacturer: Ability Dynamics</p>	NO	22 inHg	No sleeve required if used with seal-in liner	YES	INTEGRATED	YES	YES
<a href="http://rushfoot.com">rushfoot.com</a>							
 <p><b>LimbLogic</b> Manufacturer: WillowWood</p>	YES	20 inHg	Recommended LimbLogic Sleeve	YES Fresh water only	IN-LINE OR SIDE-MOUNT	YES	YES
<a href="https://www.willowwoodco.com/products-services/elevated-vacuum/limblogic/">https://www.willowwoodco.com/products-services/elevated-vacuum/limblogic/</a>							
 <p><b>Unity</b> Manufacturer: Ossur</p>	NO	15-18inHg (max 22inHg)	No sleeve but for use with Ossur Iceross Seal-In XTF	WEATHERPROOF No submersion, Salt water voids warranty	INTEGRATED	YES	YES
<a href="https://assets.ossur.com/library/35261">https://assets.ossur.com/library/35261</a>							
 <p><b>Harmony Systems</b> Manufacturer: OttoBock</p>	Harmony E2 Yes Others No	P2 - 23 inHg HD - 24 P3 - 20 inHg Triton - 20 inHg 3R60 - 21 E2 - 21 inHg	Recommended OttoBock ProFlex Plus Sleeve	E2 - WATER RESISTANT All others - YES	IN-LINE	P2, 3R60, and HD serviced by OttoBock loaner available All others YES	Triton, P3 and P4 - Yes
<a href="https://professionals.ottobockus.com/media/pdf/13338B_Harmony_Pump_Chart.pdf">https://professionals.ottobockus.com/media/pdf/13338B_Harmony_Pump_Chart.pdf</a>							
 <p><b>EchelonVac</b> Manufacturer: Endolite</p>	NO	12-17 inHg	Silcare Breathe liner	WEATHERPROOF No showering or swimming.	INTEGRATED	YES	YES
<a href="http://www.endolite.com/products/echelonvac">http://www.endolite.com/products/echelonvac</a>							

Call us toll-free at 855.450.7300 or visit us at [rushfoot.com](http://rushfoot.com)







 #ElevateYourLife

**ABILITY DYNAMICS**  
Proudly handcrafted in the U.S.A.