

Sgt. Matthew Pennington
BK, RUSH LoPro®

RUSH LoPro®

INNOVATION YOU CAN FEEL

Thanks to the unique sole plate and virtually indestructible glass composite, the RUSH LoPro® delivers the smoothest transition from heel strike to toe-off of any prosthetic foot in the industry. Patients will experience unparalleled comfort, range of motion and zero dead spots in every step.

“BEFORE RUSH®, I BROKE CARBON FEET CONSTANTLY. I BEAT MY RUSH FOOT® 24/7. IT WON’T BREAK!”

– Katie Mettner, BK, RUSH LoPro®



RUSH LoPro® EVAQ8 PROFILE INFORMATION

RUSH LoPro® EVAQ8 PROFILE INFORMATION				BUILD HEIGHT	
Weight Limit	300lbs./136kg	Heel Rise	3/8”/10mm	Size	Build Height*
Size	22cm - 29cm	Warranty	36 Months	22 - 24cm	4.5”
Foot Shell	Light or Dark			25 - 27cm	4.5”
Product Weight	649g/1.43lbs. (26cm)			28 - 29cm	4.5”

† The listing of this code should not be considered a guarantee for coverage or a payment. * Build height includes footshell

CATEGORY CHART

Weight	lbs	90 - 116	117 - 146	147 - 180	181 - 218	219 - 261	262 - 300
	kg	40 - 52	53 - 66	67 - 81	82 - 98	99 - 118	119 - 136
Low Activity		1	2	3	4	5	6
Moderate Activity		1	2	3	4	5	6
High Activity		2	3	4	5	6	7
High Impact		2-3	3-4	4-5	5-6	6-7	7-8

NOTE: Increasing category numbers indicate increased keel response.

Call us toll-free at 855.450.7300 or visit us at proteorusa.com

f t You Tube in p #HumanFirst



HumanFirst