



Erin Deegan Conner, AK, RUSH RAMPAGE LP

# INNOVATION YOU CAN FEEL

Flexibility, energy return, smooth roll through, and zero dead spots are just a few characteristics of the **RUSH RAMPAGE LP**. Packaged in a medium build-profile, the **RUSH RAMPAGE LP** is the perfect option for those with limited clearance; seeking the ultimate blend of comfort, performance, and durability.



### RUSH RAMPAGE LP PROFILE INFORMATION

Weight Limit	365lbs./165kg	Heel Rise	3/8" / 10mm	Size	Build Height**
Size	22cm - 29cm	Warranty	36 Months	22 - 24cm	4.5" / 114mm
Foot Shell	Light or Dark	†Suggested Codes	L5981, L5986	25 - 27cm	4.75" / 121mm
Product Weight*	805g			28 - 29cm	4.75" / 121mm

† The listing of this code should not be considered a guarantee for coverage or a payment.

\* Product weight based on size 26cm, Cat 4 foot module with spectra sock and foot shell.

\*\* Build height based on size 23cm, 26cm, or 29cm Cat 4 foot module with spectra sock, foot shell, and 10mm heel rise.

### CATEGORY CHART

Weight	lbs	0 - 105	106 - 140	141 - 175	176 - 210	211 - 245	246 - 280	281 - 315	316 - 350	351 - 365
	kg	0 - 48	49 - 64	65 - 79	80 - 95	96 - 111	112 - 127	128 - 143	144 - 159	160 - 166
Low Activity		1	2	3	4	5	6	7	8	9
Moderate Activity		1	2	3	4	5	6	7	8	9
High Activity		2	3	4	5	6	7	8	9	-
High Impact		3	4	5	6	7	8	9	-	-

NOTE: Increasing category numbers indicate increased keel response. FOR BILATERAL PATIENTS: Increase by one category for added stability.

Call us toll-free at 855.450.7300 or visit us at [proteorusa.com](http://proteorusa.com)

f t YouTube in p #HumanFirst



Human First