

RUSH™

RAMPAGE



Unbridled Capabilities.

The latest addition to the family, the **RUSH Foot RAMPAGE**, has arrived. We've taken the smooth roll-through, increased comfort, and all-terrain capability of its predecessors and created a new, more cosmesis-friendly option, with a lower overall product weight. The **RUSH Foot RAMPAGE** is the perfect solution for individuals seeking a blend of both cosmetic, and high-performance capabilities found only in the RUSH Foot Collection. Are you ready to **RAMPAGE?**

- Cosmesis-ready to meet your aesthetic preferences.
- Lower overall product weight.
- Smooth roll through with unparalleled energy return.
- Biometric ankle motion simulates natural gait.
- Extreme reliability for any use in any environment.
- Manufacturer recommended coding based on AOPA testing guidelines L5987 and L5986.

RUSH RAMPAGE PROFILE INFORMATION

Weight Limit	365lbs./165kg	Heel Rise	3/8" / 10mm	Size	Build Height**
Size	22cm - 29cm	Warranty	36 Months	22 - 24cm	6.625" / 168mm
Foot Shell	Light or Dark	Suggested Codes	L5987, L5986	25 - 27cm	7.125" / 181mm
Product Weight*	850g			28 - 29cm	7.25" / 184mm

† The listing of this code should not be considered a guarantee for coverage or a payment.
 * Product weight based on size 26cm, Cat 4 foot module with spectra sock and foot shell.
 **Build height based on size 23cm, 26cm, or 29cm Cat 4 foot module with spectra sock, foot shell, and 10mm heel rise.

CATEGORY CHART

Weight	lbs	0 - 105	106 - 140	141 - 175	176 - 210	211 - 245	246 - 280	281 - 315	316 - 350	351 - 365
	kg	0 - 48	49 - 64	65 - 79	80 - 95	96 - 111	112 - 127	128 - 143	144 - 159	160 - 166
Low Activity		1	2	3	4	5	6	7	8	9
Moderate Activity		1	2	3	4	5	6	7	8	9
High Activity		1	2	3	4	5	6	7	8	9
High Impact		2	3	4	5	6	7	8	9	-

NOTE: Increasing category numbers indicate increased keel response. FOR BILATERAL PATIENTS: Increase by one category for added stability.

Call us toll-free at 855.450.7300 or visit us at proteorusa.com

f t YouTube in p #HumanFirst



HumanFirst